pedalo®

Pedalo® Federbrett 150 & 180

Pedalo® Spring Board 150 & 180

130 05150

130 05180

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
The large Pedalo-Spring Board is the perfect training station for sports and therapy to improve the overall posture and movement coordination.

The platforms of 150x40 cm or 180x60 cm offer additional dynamic balance exercises from head to toe - you won't succeed comparably, on hardly any other training equipment. Ideally used for walk training (fall prevention, safe standing) - also after injuries or surgeries e. g. in the ankle, knee or hip joints. For stabilizing exercises of the entire musculoskeletal system in standing, sitting or lying positions. The unstable platform increases the number of input information (stimulus through the senses). Reaction, balance skills and especially the depth perception (proprioceptive training) is trained particularly and thereby improves the interaction between nerves and muscles. The unique spring system intelligently adapts to the body weight and also reacts to each movement of your body.

The combination of Pedalo-Textile-Powerband and Spring Board creates a great fitness station used for a modern comprehensive stabilization and coordination training on an unstable platform. The Pedalo-Textile-Powerband can be attached individually along the board and therefore allows comprehensive exercises for strength, endurance and mobility.

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
Spring Board 150 (130 05150)
 The spring board 150 is mainly suitable for exercises while standing, sitting and walking (depending on user's height).

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Spring Board 180 (130 05180)
 The area of 180x60 allows additional exercises in push-up and plank position, from the supine and prone position, as well as from the all-fours position. Perfect for functional training and therapy with adults.

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Specification:

- Multiplex birch wood 21 mm
- Surface in clear lacquer
- Integrated senso-spring system
- 4 springs fix mounted
- Spring Board 150: 21 kg
- Spring Board 180: 35 kg

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Accessories: Pedalo-Spring Element (140 06000)
 Individual increase of the spring board stability.

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Alternately stretch the hand in extension.

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Alternately stretch the leg in extension.

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Alternately stretch out the hand and leg diagonally.

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Variant:

Repeat all exercises once again but now only balance out with the knee (foot does not touch the board).

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Keep your back straight. Then push your back up as far as possible (hunchback). While doing so, bend your head towards the breastbone.

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Go into the plank position. Keep your back straight. Then move one leg alternately backwards and stabilize the position.

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Keep your back straight. Alternately extend one leg backwards and return to its initial position so that the knee is lead as close as possible to the diagonally rested elbow.

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Bend legs, lift upper body. Then move the left bent arm to the middle of the body. At the same time pull the right knee. Knee and forearm touch in the middle. Move the arm and leg back to the starting position and repeat the exercise.

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Repeat exercise with the right arm and left leg.

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From the back plank, push the pelvis horizontally upwards. Shoulder, pelvis and knee form a straight line if possible. Hold position.

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Starting from the above position: lift the leg slightly.

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Variant:

In addition, relieve one arm and stretch to the side. Then stretch in extension to the leg.

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Lift the pelvis sideways. Leg axis, pelvis and upper body form a straight line.

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Starting from the above position: spread out the leg.

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Variant:
Starting from the above position: alternately pull and stretch out the leg. Extend the arm over the head in extension.

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Two-footed position and maintain body balance.

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4. Repeat all exercises and variants in the one-leg stand.

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Two-footed position and maintain body balance. Slowly and carefully move forward over the spring board.

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2. Turn your head to one side while walking

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Move to the step position, bend the back knee and balance the position.

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Starting from the above position: jump off and perform a step change. Return softly absorbed and balance position.

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Both feet next to the spring board, jump off and return softly from above vertically on the spring board. Stabilize position and continue with low jump on the floor.

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Stabilize standing leg (leg axis), spread pendulum leg to side.

Starting from the above position: in addition do a squat with the standing leg.

Variant:

Press off the one-legged squat, jump off and return softly with the other leg.

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Stabilize the standing leg (leg axis), lead the pendulum leg up against the resistance and lift the knee. Stabilize position. Arms swing in opposite direction.

Variant:

Stretch the pendulum leg against the resistance from the front to the back until the gluteal muscle is fully tensioned.

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Stabilize the standing leg (leg axis), swing by the pendulum leg against the resistance from outside to inside of the leg (adduction).

Starting from the above position: lead the pendulum leg from the inside to the outside (abduction).

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
Stretch the arms forward against the resistance of the Textile-Powerbands and then return them to the starting position.

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Stabilize the standing leg (leg axis), stretch the pendulum leg back against the resistance until the gluteal muscles are fully tensioned. At the same time (upper body is slightly bent) move the arms upwards in the extension of the body axis.

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Pull the hands in bent position to the belly. Pull the arms back.

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Move the arms in a slightly bent posture to the breastbone. Biceps exercise.

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Move the arms in a slightly bent posture from the breastbone upwards and extend them almost completely.

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Further information about the topic 'Whole Body Training' on our website www.pedalo.de or in the theme flyer 'There is more behind it!'.

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Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

Please note:

WARNING! The device is only to be used for the purposes described in this user manual.

WARNING! This device is a very demanding one! Especially when used for the first times and/or when doing exercises in standing position, a special degree in caution, attention, circumspection and a certain degree in motor abilities is assumed.

WARNING! Place the unit only on flat level ground. Protect item for damp and wet.

This device is meant for using from persons at the age of 3 up.

Maximum weight capacity approx. 120 kg.

CE Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

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



Like Pedalo on Facebook and always be well informed. On You-Tube Pedalo-product films set you in motion.

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